

140000140918

**Solutions for Independent Living**

DLF Helpline - Telephone 0845 130 9177, Textphone 020 7432 8009

[< Back to previous page](#) |
 [Edit your report](#) |
 [Print](#) |
 [Copy to clipboard](#)

For further

information about the items listed below contact the suppliers who will be able to send you illustrated leaflets, current prices lists and other information. We recommend you look for suppliers who are members of the British Healthcare Trades Association (BHTA). Members of the BHTA sign up to a code of practice designed to eliminate 'cowboy' behavior and ensure the public can trust that members will give a good service.

Inclusion in this report does not imply any recommendation from the Disabled Living Foundation (DLF). No items have been tested by the DLF and all information is provided without legal responsibility.

Products

WHEELCHAIRS - POWERED, SCOOTERS & BUGGIES

What we include:



Liteway Scooter 3

Class 2 scooter with fold down tiller. Dismantles for storage and transportation.

Features

- three wheels
- one double front wheel
- slip-resistant foot platform
- fold down tiller
- basket
- controls
- detachable battery
- padded seat
- adjustable armrests

Accessories

- front basket
- swivel seat
- extra battery pack

Product dimensions

Function/Use

Indoor:	Standard
Outdoor:	Standard
Pavement use:	Standard

Drive

Front drive (F); rear drive (R): R

Dimensions

Dual speed option Y/N: N
 Speed (km/h): 6.4km/h

Critical dimension

Capacity: 136kg
 Length: 104cm
 Speed: 6.4km/h
 Width (cm): 56cm

Chassis

Three or four wheels: 3

Backrest

Reclining: -
 Headrest: -
 Fold down backrest: -

Accessories

Oxygen cylinder holder: -
 Legrest: -

General

Weight: 48kg
 Range: 16km
 Diameter: 18cm;
 23cm

Supplier**BHTA****Supplier's name for product****Price guide**

All Handling (Movability) Ltd
 , SW20 8DX
Tel: (020) 8542 2217
Fax: (020) 8395 4410
Email: info@movability.com



RASCAL LITEWAY SCOOTER 995

Electric Mobility Euro Ltd
 , Somerset , TA19 9DL
Tel: (01460) 258100
Fax: (01460) 258125
Email: sales@electricmobility.co.uk



Liteway Scooter 2,095-3,825

Fleinns Medicare Ltd
 , Hampshire , SO50 5PQ
Tel: (023) 8062 0060
Fax: (023) 8065 0230
Email: simon@fleinnsmedicare.co.uk



RASCAL LITEWAY SCOOTER POA

Laybrook.com
 , BS6 5QZ
Tel: (01869) 278504

LITEWAY 835

Email: sales@laybrook.com

LifeStyle Access & Mobility

, W3 9RH

Tel: (0845) 122 1486

Fax: (020) 8992 6336

Email: info@lifestyleaccess.com

LITEWAY 4

POA



PRO RIDER ROAD KING MOBILITY SCOOTER

Four-wheeled scooter which dismantles for transportation. Seat can be adjusted forwards, backwards and to a forward tilting position.

Features

- anti-tip devices
- electro-magnetic braking
- adjustable height, padded, swivel captain's seat
- width adjustable, fold-up armrests
- adjustable angle, fold-down backrest
- adjustable bar tiller
- detachable basket
- horn
- lights and indicators
- free wheel mechanism
- two batteries with direct access
- battery charge indicator

Product dimensions

Terrain

Kerb climbing: Standard

Seat Unit

Swivel seat: Standard

Seat unit

Fold-up or swing away armrests: Standard

Function/Use

Outdoor: Standard

Pavement use: Standard

Function/use

Road use: Standard

Dimensions

Speed (km/h): 12km

Length (cm): 140cm

Capacity (kg): 160kg
Width (cm): 64cm

Chassis

Three or four wheels: 4

Backrest

Fold down backrest: Standard

Supplier

BHTA

Supplier's name for product

Price guide

Pro Rider Mobility Ltd
, Leicestershire , LE2
5JG
Tel: (0116) 271 9998
Fax: (0116) 260 7733



PRO RIDER ROAD KING MOBILITY
SCOOTER

849



Rascal

Class 2 and 3 scooter. Dismantles for storage and transportation.

Features

- four wheels
- narrow base
- pneumatic tyres
- central footplate
- slip-resistant matting
- front bumper
- adjustable position tiller
- handgrips
- headlight
- central controls with one-hand forward and reverse lever
- shopping basket
- rear battery pack
- on-board charging
- rear lights
- adjustable height, padded upholstery
- captain's seat on slide base with headrest
- fold up and adjustable width armrests
- optional lightweight seat

Product dimensions

Terrain

Kerb climbing: -

Seat unit

Fold-up or swing away armrests: Standard

Seat Unit

Adjustable height: Standard
Swivel seat: Standard

Function/Use

Pavement use: Standard
Indoor: Standard
Outdoor: Standard

Dimensions

Length (cm): 125
Speed (km/h): 6.4; 9.7
Capacity (kg): 136
Dual speed option Y/N: Y
Gradient (degrees): 12
Weight of heaviest component (kg): 20
Width (cm): 54

Chassis

Three or four wheels: 4

Backrest

Fold down backrest: Standard

Accessories

Legrest: -
Oxygen cylinder holder: -

Supplier

BHTA

Supplier's name for product

Price guide

All Handling (Movability) Ltd
, SW20 8DX
Tel: (020) 8542 2217
Fax: (020) 8395 4410
Email: info@movability.com



RASCAL 388

1,995

Allardyce Healthcare Ltd
, DD4 8XD , Scotland
Tel: (01382) 773600
Fax: (01382) 774864
Email: mikes@allardycehealthcare.co.uk



RASCAL 388 SCOOTER

POA

Dolphin Stairlifts
, Dorset , BH15 2AA
Tel: (01202) 669194
Fax: (01202) 673708
Email: headoffice@dolphinstairlifts.com

RASCAL 388

1,900

Electric Mobility Euro Ltd

, Somerset , TA19 9DL

Tel: (01460) 258100**Fax:** (01460) 258125**Email:** sales@electricmobility.co.uk

Rascal

2,875

Euromobility Ltd

, Kent , ME16 0PU

Tel: (01622) 755744**Fax:** (01622) 766233**Email:** info@euromedical.co.uk

RASCAL 388

899

Factory Outlet Scooters

, Lincolnshire , NG31 6EJ

Email: protectmobility@aol.comELECTRIC MOBILITY RASCAL
388

749

Fleinns Medicare Ltd

, Hampshire , SO50 5PQ

Tel: (023) 8062 0060**Fax:** (023) 8065 0230**Email:** simon@fleinnsmedicare.co.uk

RASCAL 388

POA

Laybrook.com

, BS6 5QZ

Tel: (01869) 278504**Email:** sales@laybrook.com

RASCAL 388

1,045

LifeStyle Access & Mobility

, W3 9RH

Tel: (0845) 122 1486**Fax:** (020) 8992 6336**Email:** info@lifestyleaccess.comELECTRIC MOBILITY RASCAL
388

1,995

**Revo 3**

Class 2 scooter. Dismantles for storage and transportation.

Features

- three or four wheels
- slip-resistant foot platform
- rear wheel drive
- adjustable fold down tiller
- lights
- adjustable height swivel seat
- fold up armrests
- backrest
- on-board accessory bracket
- free-wheel lever
- charger point and circuit breaker

Accessories

- lap belt

- stick and crutch holder
- cup holder
- basket
- frame holder
- stick clip attachment

Supplier**BHTA****Discount Mobility Plus Ltd**

, Rutland , LE15 6AY

Tel: (0808) 178 4151**Fax:** (01572) 801075**Email:**

sales@discountmobilityplus.co.uk

Easy Mobility Ltd

, Surrey , KT15 2AF

Tel: (0800) 626829**Email:** info@easymobility.co.uk**Factory Outlet Scooters**

, Lincolnshire , NG31 6EJ

Email: protectmobility@aol.com**Laybrook.com**

, BS6 5QZ

Tel: (01869) 278504**Email:** sales@laybrook.com**LifeStyle Access & Mobility**

, W3 9RH

Tel: (0845) 122 1486**Fax:** (020) 8992 6336**Email:** info@lifestyleaccess.com**Morecare Mobility & Healthcare Specialists**

, Wiltshire , SN10 2AG

Tel: (01380) 727555**Fax:** (01380) 736696**Email:**

info@morecarehealthcare.co.uk

Pride Mobility Products Ltd

, Oxfordshire , OX26 4UL

Tel: (01869) 324600**Fax:** (01869) 323070**Simplyhealth**

, Hampshire , SP10 5AR

Tel: (01264) 350022**Fax:** (01264) 350033**Email:**

info@simplyhealthstore.co.uk





ROMA PARIS GK9

Four-wheeled scooter which dismantles for storage and transportation.

Features

- height adjustable, full swivel seat which locks in four positions
- height and width adjustable arms
- fold up padded armrests
- fold down backrest
- multi-adjustable tiller angle
- horn and basket

Product dimensions

Seat unit

Fold-up or swing away armrests: Standard

Seat Unit

Swivel seat: Standard

Function/Use

Outdoor: Standard

Pavement use: Standard

Function/use

Road use: -

Dimensions

Weight of heaviest component (kg): 16kg

Width (cm): 56cm

Speed (km/h): 6

Gradient (degrees): 8

Length (cm): 105cm

Capacity (kg): 136kg

Dual speed option Y/N: N

Chassis

Three or four wheels: 4

Backrest

Fold down backrest: Standard

Supplier

BHTA

Supplier's name for product

Price guide

All Handling (Movability) Ltd
, SW20 8DX
Tel: (020) 8542 2217
Fax: (020) 8395 4410
Email: info@movability.com



ROMA SHOPRIDER PARIS GK9 POA

DCS Joncare
, Oxon , OX14 3RY
Tel: (01235) 523353
Fax: (01235) 531019
Email: info@dcsjoncare.com



ROMA PARIS GK9 POA

RDK Mobility
, Manchester , M26 1NR
Tel: (0161) 724 9580
Fax: (0161) 724 4453
Email: sales@rdkmobility.co.uk



SHOPRIDER PARIS 21040810 795-945

Roma Medical Aids Ltd
, CF31 3TB
Tel: (01656) 674488
Fax: (01656) 674499
Email: sales@romamedical.co.uk



ROMA PARIS GK9 POA



STERLING PEARL

Class 2 scooter. Dismantles for storage and transportation.

Features

- rear wheel drive
- four wheels
- gem-lock mechanism
- slip-resistant foot platform
- bumpers
- adjustable seat height
- desk style armrests
- padded backrest
- controls
- front basket
- detachable batteries

Product dimensions

Seat Unit

Adjustable height: Standard

Seat unit

Fold-up or swing away armrests: Standard

Seat Unit

Swivel seat: Standard

Function/Use

Indoor: Option
Outdoor: Standard
Pavement use: Standard

Dimensions

Speed (km/h): 6
Length (cm): 103
Gradient (degrees): 8
Capacity (kg): 136
Dual speed option Y/N: N
Turning circle (cm): 117
Weight of heaviest component (kg): 16
Width (cm): 56

Chassis

Three or four wheels: 4

Accessories

Oxygen cylinder holder: -
Legrest: -

Supplier

BHTA

Supplier's name for product

Price guide

Care Necessities

, Staffordshire , WS15 4PT
Tel: (01543) 493235
Fax: (01543) 493485
Email: info@carenecessities.co.uk



STERLING PEARL

POA

Factory Outlet Scooters

, Lincolnshire , NG31 6EJ
Email: protectmobility@aol.com



SUNRISE MEDICAL STERLING PEARL

735

Fleinns Medicare Ltd

, Hampshire , SO50 5PQ
Tel: (023) 8062 0060
Fax: (023) 8065 0230
Email: simon@fleinnsmedicare.co.uk



STERLING PEARL SCOOTER

POA

Helpful Hand

, Cheshire , SK11 8DG
Tel: (01625) 424438
Fax: (01625) 617857
Email: info@thehelpfulhand.co.uk



STERLING PEARL SCOOTER

POA

John Bell & Croyden

, W1U 2AU

Tel: (020) 7935 5555**Fax:** (020) 7935 9605**Email:** jbc@johnbellcroyden.co.uk

STERLING PEARL

POA

LifeStyle Access & Mobility

, W3 9RH

Tel: (0845) 122 1486**Fax:** (020) 8992 6336**Email:** info@lifestyleaccess.com

STERLING PEARL

POA

Mobility Zone Ltd

, LU7 4UH

Tel: (01525) 382179**Fax:** (01525) 378347**Email:**

sales@mobilitymegastore.co.uk

SUNRISE MEDICAL MOBILITY
PEARL SCOOTER

725

**Nightingale Medical & Mobility
Products**

, Berkshire , SL1 7PE

Tel: (01628) 668660**Fax:** (01628) 663584

SMALL SCOOTER 6184

1,375

Sunrise Medical Ltd

, West Midlands , DY8 4PS

Tel: (01384) 446688**Fax:** (01384) 446699**Email:** sunmail@sunmed.co.uk

STERLING PEARL SCOOTER

POA

The Force Ten Co Ltd

, Surrey , KT23 4RQ

Tel: (01372) 450887**Fax:** (01372) 450771**Email:** sales@forcetenco.co.uk

STERLING PEARL

995

Upson Mobility

, Essex , CM15 0UH

Tel: (01277) 352000**Fax:** (01277) 354761**Email:** info@upsonmobility.co.uk

STERLING PEARL

995

**STERLING SAPPHIRE LS SCOOTER**

Compact Class 2 scooter.

Features

- rear wheel drive
- slip-resistant foot platform
- bumpers
- four wheels
- lights
- front basket

- padded swivel seat
- adjustable fold up armrests

Product dimensions**Seat Unit**

Adjustable height: Standard

Seat unit

Fold-up or swing away armrests: Standard

Seat Unit

Swivel seat: Standard

Function/Use

Indoor: Standard

Outdoor: Standard

Dimensions

Speed (km/h): 6
 Length (cm): 122
 Gradient (degrees): 11;
 14
 Capacity (kg): 137
 Dual speed option Y/N: N
 Turning circle (cm): 140
 Width (cm): 53

Chassis

Three or four wheels: 4

Supplier**BHTA****Supplier's name for product****Price guide****Care Necessities**

, Staffordshire , WS15 4PT

Tel: (01543) 493235

Fax: (01543) 493485

Email: info@carenecessities.co.uk



STERLING SAPPHIRE LS
SCOOTER

POA

Factory Outlet Scooters

, Lincolnshire , NG31 6EJ

Email: protectmobility@aol.com



SUNRISE MEDICAL SAPPHIRE LS 790

Mobility Zone Ltd

, LU7 4UH

Tel: (01525) 382179

Fax: (01525) 378347

Email:

sales@mobilitymegastore.co.uk

SUNRISE MEDICAL STERLING
SAPPHIRE LS SCOOTER

POA

Morecare Mobility & Healthcare Specialists

, Wiltshire , SN10 2AG

Tel: (01380) 727555**Fax:** (01380) 736696**Email:**

info@morecarehealthcare.co.uk

STERLING SAPPHIRE LS 4
WHEELED SCOOTER SAPPLS

POA

Nightingale Medical & Mobility Products

, Berkshire , SL1 7PE

Tel: (01628) 668660**Fax:** (01628) 663584

SMALL SCOOTER 6163

1,550

Sunrise Medical Ltd

, West Midlands , DY8 4PS

Tel: (01384) 446688**Fax:** (01384) 446699**Email:** sunmail@sunmed.co.ukSTERLING SAPPHIRE LS
SCOOTER

POA

The Force Ten Co Ltd

, Surrey , KT23 4RQ

Tel: (01372) 450887**Fax:** (01372) 450771**Email:** sales@forcetenco.co.ukSTERLING SAPPHIRE LS
SCOOTER

POA

Upson Mobility

, Essex , CM15 0UH

Tel: (01277) 352000**Fax:** (01277) 354761**Email:** info@upsonmobility.co.uk

STERLING SAPPHIRE LS SAPHY8 1,195

**Sterling Sapphire Scooter**

Class 2 outdoor scooter.

Features

- three or four wheels
- pneumatic tyres
- rear wheel drive
- batteries
- slip-resistant foot platform
- front and rear lights
- indicators
- adjustable tiller
- charging point
- steering handles with thumb pressure accelerator and brake levers
- adjustable speed control
- automatic release braking
- adjustable height and width rotating padded sliding captain seat
- headrest
- fold up desk style armrests
- optional basket

Product dimensions

Seat Unit

Adjustable height:	Standard
Swivel seat:	Standard

Function/use

Road use: -

Function/Use

Indoor:	Standard
Outdoor:	Standard
Pavement use:	Standard

Dimensions

Speed (km/h):	6
Length (cm):	122
Capacity (kg):	136
Dual speed option Y/N:	N
Gradient (degrees):	14
Turning circle (cm):	74; 111
Weight of heaviest component (kg):	21.5; 27
Width (cm):	53

Chassis

Three or four wheels: 3; 4

Supplier

BHTA

Supplier's name for product

Price guide

All Handling (Movability) Ltd
 , SW20 8DX
Tel: (020) 8542 2217
Fax: (020) 8395 4410
Email: info@movability.com



SAPHIRE 2 SCOOTER

1,100

Care Necessities
 , Staffordshire , WS15 4PT
Tel: (01543) 493235
Fax: (01543) 493485
Email: info@carenecessities.co.uk



Sterling Sapphire Scooter

POA

LifeStyle Access & Mobility
 , W3 9RH
Tel: (0845) 122 1486
Fax: (020) 8992 6336
Email: info@lifestyleaccess.com

STERLING SAPPHIRE SCOOTER

POA

Mobility Zone Ltd
 , LU7 4UH
Tel: (01525) 382179
Fax: (01525) 378347

SUNRISE MEDICAL STERLING

995

Email:

sales@mobilitymegastore.co.uk

Powerchairs Mobility Consultants

, B74 2HE

Tel: (0121) 353 7975**Fax:** (0121) 353 7978**Email:** sales@power-chairs.com

STERLING SAPPHIRE

POA

Prestige Healthcare (London) Ltd

, Hertfordshire , EN5 5BX

Tel: (0870) 745 2300**Fax:** (020) 8441 9123**Email:** lornastump@yahoo.co.ukSTERLING SAPPHIRE 4-WHEEL
SCOOTER SAP4

1,995

Sunrise Medical Ltd

, West Midlands , DY8 4PS

Tel: (01384) 446688**Fax:** (01384) 446699**Email:** sunmail@sunmed.co.uk

STERLING SAPPHIRE SCOOTER

POA

The Force Ten Co Ltd

, Surrey , KT23 4RQ

Tel: (01372) 450887**Fax:** (01372) 450771**Email:** sales@forcetenco.co.uk

STERLING SAPPHIRE SCOOTER

POA

Upson Mobility

, Essex , CM15 0UH

Tel: (01277) 352000**Fax:** (01277) 354761**Email:** info@upsonmobility.co.uk

STERLING SAPPHIRE 3 WHEEL

POA

**STRIDER MIDI SCOOTER**

Class 2 scooter. Range of colours available. Dismantles for storage and transportation.

Features

- three or four wheels
- slip-resistant foot platform
- rear wheel drive
- lights
- sliding adjustable swivel seat
- adjustable width and height armrests
- push-pull tiller controls
- mirror
- front and rear bumpers
- integral lifting handle
- front and rear lights
- automatic brake
- manual hand brake on MD4

Accessories

- oxygen holder
- crutch and cane holder
- rear basket
- captain seat

Product dimensions

Terrain

Kerb climbing: Standard

Seat Unit

Adjustable height: Standard

Seat unit

Fold-up or swing away armrests: Standard

Seat Unit

Swivel seat: Standard

Function/Use

Indoor: Standard

Outdoor: Standard

Pavement use: Standard

Dimensions

Speed (km/h): 6

Weight of heaviest component (kg): 24

Width (cm): 56

Turning circle (cm): -

Length (cm): 124

Gradient (degrees): 8-10

Capacity (kg): 150

Dual speed option Y/N: N

Chassis

Three or four wheels: 3; 4

Backrest

Fold down backrest: -

Accessories

Oxygen cylinder holder: Option

Legrest: Standard

Supplier

BHTA

Supplier's name for product

Price guide

All Handling (Movability) Ltd
 , SW20 8DX
Tel: (020) 8542 2217
Fax: (020) 8395 4410
Email: info@movability.com



DMA STRIDER MIDI SCOOTER POA

Factory Outlet Scooters
 , Lincolnshire , NG31 6EJ
Email: protectmobility@aol.com



DMA DAYS MEDICAL STRIDER
 MIDI 3 995

Helpful Hand
 , Cheshire , SK11 8DG
Tel: (01625) 424438
Fax: (01625) 617857
Email: info@thehelpfulhand.co.uk



DMA MIDI 3 POA

Nightingale Medical & Mobility Products
 , Berkshire , SL1 7PE
Tel: (01628) 668660
Fax: (01628) 663584

MEDIUM SCOOTER 632-4 1,845

Rose Health Care
 , Suffolk , IP1 3EP
Tel: (01473) 258508
Fax: (01473) 258508
Email: rosehealthcare@btconnect.com



MIDI R86 1,250

Class 2 scooters with fold down tiller

What we include: Included in this section are compact Class 2 scooters with tillers which fold down to assist storage or transportation. They may be three or four-wheeled and most models are not intended to climb kerbs. This range includes lighter weight scooters, intended for occasional or semi-permanent use. They will not usually have lights fitted, but may include accessories such as shopping baskets.



DMA STRIDER MINI SCOOTER

Class 2 indoor/outdoor scooter. Dismantles for storage and transportation.

Features

- three or four wheels
- rear wheel drive: fold down tiller
- slip-resistant foot platform
- anti-tippers
- front basket
- padded adjustable height swivel seat
- adjustable width and angle padded desk style armrests
- detachable battery pack
- front and rear carry handles
- rear bumper
- mirror
- hand controls
- on-board battery charger
- free-wheel lever

Accessories

- Delta tiller kit

- crutch and stick holder
- rear basket
- raincoat for scooter and user
- saddle bag
- oxygen bottle holder
- walking frame holder
- long range battery pack

Product dimensions

Critical dimension

Seat Depth:	43cm
Seat Height:	39-47cm
Seat Width:	44cm
Capacity:	130kg
Length:	110cm
Height:	91cm
Speed:	6km
Width (cm):	49cm

General

Range:	21km
Seat Height:	49-57cm

Supplier	BHTA	Supplier's name for product	Price guide
All Handling (Movability) Ltd , SW20 8DX Tel: (020) 8542 2217 Fax: (020) 8395 4410 Email: info@movability.com		DMA STRIDER MINI SCOOTER	850
Factory Outlet Scooters , Lincolnshire , NG31 6EJ Email: protectmobility@aol.com		DMA DAYS MEDICAL STRIDER MINI 3	849
GBL Wheelchair Services Ltd , Middlesex , TW13 6AR Tel: (020) 8867 2461 Fax: (020) 8867 2475 Email: info@gbwheelchairs.com		STRIDER MINI 4 SCOOTER	875
Helpful Hand , Cheshire , SK11 8DG Tel: (01625) 424438 Fax: (01625) 617857 Email: info@thehelpfulhand.co.uk		DMA MINI 3	POA
LifeStyle Access & Mobility , W3 9RH Tel: (0845) 122 1486 Fax: (020) 8992 6336 Email: info@lifestyleaccess.com		DAYS MEDICAL STRIDER MINI 3 AND 4 WHEELER	1,350
Nightingale Medical & Mobility Products , Berkshire , SL1 7PE Tel: (01628) 668660 Fax: (01628) 663584		SMALL SCOOTER 6204	1,395

Simplyhealth

, Hampshire , SP10 5AR

Tel: (01264) 350022**Fax:** (01264) 350033**Email:** info@simplyhealthstore.co.uk

MINI SCOOTER

POA

Yorkshire Care Equipment

, West Yorkshire , LS19 6DY

Tel: (0113) 250 3250**Fax:** (0113) 250 7433**Email:** sales@yorkshirecare.com

STRIDER MINI 4

POA

Subscriber Telephone: 020 7432 8003

[Contact Information](#) | [User policy](#) | [Privacy policy](#)

Registered office: Disabled Living Foundation, 380-384 Harrow Road, London, W9 2HU | Tel: 020 7289 6111

Disabled Living Foundation is a Company Limited by Guarantee Registered in England Number 1837993 |
Registered Charity No. 290069.

Copyright 2009 Disabled Living Foundation. All rights reserved.